APPS PATIENT POSITIONING SYSTEMS

The Agency for Healthcare Research and Quality estimates that over 700,000 hospitalized patients and nearly half of the 1.6 million nursing home residents in the United States fall each year.

The PPS LIFT® is a fall-recovery device that effortlessly lifts patients from the floor to bed or stretcher height in a supine position. Patient comfort is maximized while the risk of lifting related injury to both the patient and caregiver is minimized. The PPS LIFT is a critical component for fall prevention and emergency preparedness.



INTRODUCING THE



Promotes Patient Comfort & Dignity

- Supine lifting feature prevents patient from incurring further injuries after a fall
- Single inflation point that rapidly inflates the device WITHOUT needing to disconnect the air hose
- Exceedingly stable
- ▶ Weight capacity of 1200 lbs.








Improves Staff Efficiency & Safety

- Air-assisted technology reduces the number of caregivers required to lift a patient when responding to a fall.
- Superior to the competition
- ▶ Using the PPS LIFT in conjunction with the PPS Glide lateral air transfer system provides your caregivers with a safe no-lift solution. After the individual has been lifted using the PPS LIFT, the patient can then be laterally transfer them to a bed or stretcher.





BENEFITS OF THE **PPSLift** AIR PATIENT LIFT:

- System flexibility Allows for use in spaces where floor or ceiling lifts cannot be utilized
- ► Teflon® coated bottom Makes it easier to transfer the individual
- Radiolucent and MRI compatible –
 For artifact-free imaging
- Available in 39" (1200 lb weight limit)



PPSLift

CONTACT US







www.PPSproducts.com